



Suction Cup Accessory
5200-50-BX

Doc Johnson's® Main Squeeze strokers incorporate a hard-case construction that puts the power in the user's hands, with a squeeze plate designed to vary the tightness of each stroke, and a suction-adjusting end cap that twists to multiple levels of intensity. Each Main Squeeze™ stroker features an ULTRASKYN™ sleeve with a uniquely textured interior, for ultra-pleasurable stimulation. Main Squeeze's sleek, durable case and screw-on top keep it clean, discreet, and convenient for travel; and Main Squeeze™ strokers also fully disassemble for easy cleanup. All items from Main Squeeze™ include detailed instructions on usage and care.

In addition to Main Squeeze's selection of hard-case strokers, The Main Squeeze™ line also offers a range of curated accessories designed to enhance the Main Squeeze™ experience.

Featured Products

Suction Cup Accessory - 5200-50-BX

Attaching securely to any of Doc Johnson's Main Squeeze™ full-sized strokers, the Suction Cup Accessory powerfully adheres to any smooth, non-porous surface. The Main Squeeze Suction Cup Accessory is designed for hands-free use of any full-sized Main Squeeze™ stroker and is ideal for shower/bath use. The adjustable arm swivels and bends for custom positioning, and locks securely into place. Includes complete instructions for use.

Warming Accessory - 5204-01-BX

The Main Squeeze™ Warming Accessory heats up any Main Squeeze's ULTRASKYN™ interior to a body-safe 98°F (37°C) in just 8 minutes. The Warming Accessory was created exclusively for use with Main Squeeze™ strokers, approximating body temperature to feel even more like the real thing, and featuring a slim design that won't alter the shape, texture, or tightness of the product. Connect the attached USB cord to any free outlet before pushing the sleek shape deep into a stroker's length. As with any warming product, users should test heat level before using.

The Virgin - Insert Replacement - 5202-04-CD

This ULTRASKYN™ insert is a replacement sleeve for Doc Johnson's 'Main Squeeze™ - The Virgin' hard case masturbator. This sleeve features the same extra-tight interior and realistic breaking hymen that comes included with the full Virgin Main Squeeze™ product, and can also be used with any other Main Squeeze™ hard-case exteriors.



Warming Accessory
5204-01-BX



The Virgin Insert Replacement
5202-04-CD

Care and Maintenance

To clean Suction Cup Accessory, Warming Accessory, and The Virgin Insert Replacement, spray with Doc Johnson Toy Cleaner or wash with mild soap and warm water. Dust The Virgin Insert Replacement with Doc Johnson Refresh Powder before storing. Allow all items to dry completely before storing in a cool, dry place.

Dimensions

	Total Length	Width/Diameter	Circumference
Suction Cup Accessory	6.5" (16.51 cm)	3.5" (7.62 cm)	
Warming Accessory	7.5" (19.05 cm)		
Virgin Insert Replacement	7.5" (19.05 cm)	3" (7.62 cm)	9.5" (24.13 cm)

CONTACT INFO: